

67kDa Laminin Receptor Rabbit mAb


ABZ2811

Product Overview

Name	67kDa Laminin Receptor Rabbit mAb
Accession(Primary)	P08865
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 33 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:200-1:2000; IP-1:10-1:100

Target information (P08865)

Gene Name	RPSA
About Gene Symbol	RPSA (Ribosomal Protein SA) is an extracellular matrix laminin protein involved in translational regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of 67kDa Laminin Receptor in K562, C6, 3T3, HeLa lysates using 67kDa Laminin Receptor antibody.