

Neurofilament light polypeptide/NF-L Rabbit mAb


ABZ2812

Product Overview

Name	Neurofilament light polypeptide/NF-L Rabbit mAb
Accession(Primary)	P08551
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 62 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; ICC/IF-1:20-1:50; IP-1:20-1:50

Target information (P08551)

Gene Name	Nefl
About Gene Symbol	



Western blot analysis of 68kDa Neurofilament in mouse brain lysates using Neurofilament light polypeptide antibody.