

AAMP Rabbit mAb


ABZ2813

Product Overview


Name	AAMP Rabbit mAb
Accession(Primary)	Q13685
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:10-1:100

Target information (Q13685)

Gene Name	AAMP
About Gene Symbol	AAMP (Angio Associated Migratory Cell Protein) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of AAMP in mouse heart lysates using AAMP antibody.



Western blot analysis of AAMP in C6, HeLa lysates using AAMP antibody.