

Activin A Receptor Type IB/ALK-4 Rabbit mAb

ABZ2817


Product Overview

Name	Activin A Receptor Type IB/ALK-4 Rabbit mAb
Accession(Primary)	P36896
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 57 kDa; Observed MW: 57 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:20-1:100


Target information (P36896)

Gene Name ACVR1B

About Gene Symbol



Western blot analysis of Activin A Receptor Type IB in mouse brain, rat brain lysates using Activin A Receptor Type IB antibody.



Western blot analysis of Activin A Receptor Type IB in K562 lysates using Activin A Receptor Type IB antibody.