

Activin A Receptor Type IB/ALK-4 Rabbit mAb


ABZ2818

Product Overview


Name	Activin A Receptor Type IB/ALK-4 Rabbit mAb
Accession(Primary)	P36896
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 57 kDa; Observed MW: 57 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000; ICC/IF-1:100-1:200

Target information (P36896)


Gene Name	ACVR1B
About Gene Symbol	



Western blot analysis of Activin A Receptor Type IB in mouse brain lysates using Activin A Receptor Type IB antibody.



Western blot analysis of Activin A Receptor Type IB in CHO-K1, rat Brain lysates using Activin A Receptor Type IB antibody.



Immunocytochemistry analysis of Activin A
Receptor Type IB (green) in hela using Activin A
Receptor Type IB antibody, and DAPI(blue)