


ADH5 Rabbit mAb

ABZ2819

Product Overview

Name	ADH5 Rabbit mAb
Accession(Primary)	P11766
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:1000-1:5000

Target information (P11766)

Gene Name ADH5**About Gene Symbol**

Western blot analysis of ADH5 in MCF-7, HepG2, K562 lysates using ADH5 antibody.

