

Adrenodoxin Rabbit mAb

ABZ2820


Product Overview

Name	Adrenodoxin Rabbit mAb
Accession(Primary)	P10109
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 19 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:100; ICC/IF-1:100


Target information (P10109)

Gene Name FDX1

About Gene Symbol



Western blot analysis of Adrenodoxin in 3T3 lysates using Adrenodoxin antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Adrenodoxin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.