

AF10 Rabbit mAb


ABZ2821

Product Overview

Name	AF10 Rabbit mAb
Accession(Primary)	P55197
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 113 kDa; Observed MW: 140 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; FC-1:200-1:500

Target information (P55197)

Gene Name	MLLT10
About Gene Symbol	



Western blot analysis of AF10 in HeLa, K562
lysates using AF10 antibody.