

AGXT Rabbit mAb

ABZ2822


Product Overview

Name	AGXT Rabbit mAb
Accession(Primary)	P21549
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P21549)

Gene Name

AGXT

About Gene Symbol

Western blot analysis of AGXT in mouse liver, rat liver lysates using AGXT antibody.

