

AHA1 Rabbit mAb

ABZ2823


Product Overview

Name	AHA1 Rabbit mAb
Accession(Primary)	O95433
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 38 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:20-1:50


Target information (O95433)

Gene Name AHSA1

About Gene Symbol



Western blot analysis of AHA1 in rat brain, mouse skin lysates using AHA1 antibody.



Western blot analysis of AHA1 in HeLa, HepG2
lysates using AHA1 antibody.