

ALAS2/ASB Rabbit mAb

ABZ2824

Product Overview

Name	ALAS2/ASB Rabbit mAb
Accession(Primary)	P22557
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 65 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (P22557)

Gene Name ALAS2

About Gene Symbol



