

alpha 1a Adrenergic Receptor Rabbit mAb


ABZ2826

Product Overview

Name	alpha 1a Adrenergic Receptor Rabbit mAb
Accession(Primary)	P35348
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 51 kDa; Observed MW: 51 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:200; ICC/IF-1:200-1:500; FC-1:100-1:1000

Target information (P35348)

Gene Name	ADRA1A
About Gene Symbol	



Western blot analysis of alpha 1a Adrenergic Receptor in rat Brain, C6 lysates using alpha 1a Adrenergic Receptor antibody.