

alpha Actinin 4 Rabbit mAb


ABZ2827

Product Overview


Name	alpha Actinin 4 Rabbit mAb
Accession(Primary)	O43707
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 105 kDa; Observed MW: 105 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:20-1:50; IP-1:20-1:50

Target information (O43707)

Gene Name	ACTN4
About Gene Symbol	ACTN4 (Actinin Alpha 4) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of alpha Actinin 4 in K562, rat Brain, C6, 3T3, Hela lysates using alpha Actinin 4 antibody.



Immunocytochemistry analysis of alpha Actinin 4 (green) in Hela using alpha Actinin 4 antibody, and DAPI (blue)