

alpha Sarcoglycan Rabbit mAb


ABZ2828

Product Overview

Name	alpha Sarcoglycan Rabbit mAb
Accession(Primary)	Q16586
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000; IP-1:20-1:50

Target information (Q16586)

Gene Name	SGCA
About Gene Symbol	



Western blot analysis of alpha Sarcoglycan in mouse heart lysates using alpha Sarcoglycan antibody.