

ALS2CR8 Rabbit mAb

ABZ2829

Product Overview

Name	ALS2CR8 Rabbit mAb
Accession(Primary)	Q8N187
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 81 kDa; Observed MW: 81 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:50-1:100

Target information (Q8N187)

Gene Name CARF**About Gene Symbol**

