


Angiopoietin 2/ANG2 Rabbit mAb

ABZ2830

Product Overview

Name	Angiopoietin 2/ANG2 Rabbit mAb
Accession(Primary)	O15123
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 57 kDa; Observed MW: 57 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (O15123)

Gene Name ANGPT2**About Gene Symbol**

Western blot analysis of Angiopoietin 2 in rat Brain lysates using Angiopoietin 2 antibody.

