

Angiotensin Converting Enzyme 1 Rabbit mAb

ABZ2831


Product Overview

Name	Angiotensin Converting Enzyme 1 Rabbit mAb
Accession(Primary)	P09470
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 151 kDa; Observed MW: 180 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:500


Target information (P09470)

Gene Name Ace

About Gene Symbol



Western blot analysis of Angiotensin Converting Enzyme 1 in rat liver lysates using Angiotensin Converting Enzyme 1 antibody.



Western blot analysis of Angiotensin Converting Enzyme 1 in mouse liver, mouse heart lysates using Angiotensin Converting Enzyme 1 antibody.