

Angiotensinogen Rabbit mAb


ABZ2832

Product Overview

Name	Angiotensinogen Rabbit mAb
Accession(Primary)	P11859
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; IP-1:20-1:50

Target information (P11859)

Gene Name	Agt
About Gene Symbol	



Western blot analysis of Angiotensinogen in mouse liver lysates using Angiotensinogen antibody.