

Annexin 2 Rabbit mAb


ABZ2833

Product Overview


Name	Annexin 2 Rabbit mAb
Accession(Primary)	P07355
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 39 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:10-1:100

Target information (P07355)

Gene Name	ANXA2
About Gene Symbol	



Western blot analysis of Annexin A2 in C6, 3T3, HeLa lysates using Annexin A2 antibody.



Immunocytochemistry analysis of Annexin A2 (green) in HeLa using Annexin A2 antibody, and DAPI (blue)