

Transcription factor AP-2-alpha Rabbit mAb


ABZ2835

Product Overview

Name	Transcription factor AP-2-alpha Rabbit mAb
Accession(Primary)	P05549
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:20-1:50; FC-1:20-1:50; IP-1:20

Target information (P05549)

Gene Name	TFAP2A
About Gene Symbol	TFAP2A (Transcription Factor AP-2 Alpha) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of AP2alpha in C6 lysates
using AP2 alpha antibody.