

ApoER2 Rabbit mAb


ABZ2836

Product Overview


Name	ApoER2 Rabbit mAb
Accession(Primary)	Q14114
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 106 kDa; Observed MW: 130,106 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (Q14114)

Gene Name	LRP8
About Gene Symbol	LRP8 (LDL Receptor Related Protein 8) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of ApoER2 in Human placenta lysates using ApoER2 antibody.



Western blot analysis of ApoER2 in CHO-K1, C6 lysates using ApoER2 antibody.