


ARHGAP25 Rabbit mAb

ABZ2838

Product Overview

Name	ARHGAP25 Rabbit mAb
Accession(Primary)	P42331
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 73 kDa; Observed MW: 73 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (P42331)

Gene Name ARHGAP25**About Gene Symbol**

Western blot analysis of ARHGAP25 in rat brain, Raoms, Jurkat lysates using ARHGAP25 antibody.

