

ERAP1 Rabbit mAb

ABZ2841


Product Overview

Name	ERAP1 Rabbit mAb
Accession(Primary)	Q9NZ08
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 107 kDa; Observed MW: 107 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000


Target information (Q9NZ08)

Gene Name ERAP1

About Gene Symbol



Western blot analysis of ARTS1/ERAP1 in HeLa, C6, Jurkat lysates using ERAP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using ARTS1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.