


AS160 Rabbit mAb

ABZ2842

Product Overview

Name	AS160 Rabbit mAb
Accession(Primary)	O60343
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 147 kDa; Observed MW: 160 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (O60343)

Gene Name TBC1D4**About Gene Symbol**

Western blot analysis of AS160 in 293 lysates using AS160 antibody.

