


ATP1B1 Rabbit mAb


ABZ2843

Product Overview


Name	ATP1B1 Rabbit mAb
Accession(Primary)	P05026
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 45-55 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (P05026)

Gene Name ATP1B1**About Gene Symbol**



Western blot analysis of ATP1B1 in rat Brain lysates using ATP1B1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using ATP1B1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.