

## ATP1B3 Rabbit mAb


ABZ2844

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | ATP1B3 Rabbit mAb                          |
| <b>Accession(Primary)</b>   | P54709                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 32 kDa; Observed MW: 40 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:50-1:100         |

### Target information (P54709)

|                          |                                                                                                                                                                                                                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | ATP1B3                                                                                                                                                                                                                                                                                   |
| <b>About Gene Symbol</b> | ATP1B3 (ATPase Na <sup>+</sup> /K <sup>+</sup> Transporting Subunit Beta 3) is a cell surface marker involved in adaptive immunity. It is an important research target in immunology, flow cytometry, and cancer immunotherapy and is widely studied to elucidate cellular mechanisms an |



Western blot analysis of ATP1B3 in C6, 3T3  
lysates using ATP1B3 antibody.