


ATP5G Rabbit mAb

ABZ2845


Product Overview

Name	ATP5G Rabbit mAb
Accession(Primary)	P05496
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (P05496)

Gene Name ATP5MC1**About Gene Symbol**

Western blot analysis of ATP5G1/G2/G3 in Jurkat, C6 lysates using ATP5G antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using ATP5G1/G2/G3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.