

AZI2/NAP1 Rabbit mAb

ABZ2848


Product Overview

Name	AZI2/NAP1 Rabbit mAb
Accession(Primary)	Q9H6S1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100

Target information (Q9H6S1)

Gene Name AZI2

About Gene Symbol



Western blot analysis of AZI2 in K562, rat Brain, C6, HeLa lysates using AZI2 antibody.

