


BAT5 Rabbit mAb

ABZ2850


Product Overview

Name	BAT5 Rabbit mAb
Accession(Primary)	O95870
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 63 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (O95870)

Gene Name ABHD16A**About Gene Symbol**

Western blot analysis of BAT5 in rat muscle
lysates using BAT5 antibody.



Western blot analysis of BAT5 in Hela, Jurkat
lysates using BAT5 antibody.