

beta IV Tubulin Rabbit mAb


ABZ2851

Product Overview


Name	beta IV Tubulin Rabbit mAb
Accession(Primary)	P04350
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; ICC/IF-1:200-1:500; FC-1:50-1:100

Target information (P04350)


Gene Name	TUBB4A
About Gene Symbol	




Western blot analysis of beta IV Tubulin in mouse brain lysates using beta IV Tubulin antibody.



Western blot analysis of beta IV Tubulin in K562, rat Brain, C6, HeLa lysates using beta IV Tubulin antibody.



Immunocytochemistry analysis of beta IV Tubulin(green) in Hela using beta IV Tubulin antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human brain using beta IV Tubulin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.