

BIN1 Rabbit mAb


ABZ2853

Product Overview


Name	BIN1 Rabbit mAb
Accession(Primary)	O00499
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 65 kDa; Observed MW: 45-80 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; ICC/IF-1:50-1:100; FC-1:50-1:100

Target information (O00499)


Gene Name	BIN1
About Gene Symbol	BIN1 (Bridging Integrator 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Immunocytochemistry analysis of BIN1 (green) in U87-MG using BIN1 antibody, and DAPI (blue).



Western blot analysis of BIN1 in rat brain, mouse muscle lysates using BIN1 antibody.



Western blot analysis of BIN1 in rat Brain lysates using BIN1 antibody



Immunohistochemistry analysis of paraffin-embedded Human tonsil using BIN1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.