

BIN1 Rabbit mAb


ABZ2854

Product Overview

Name	BIN1 Rabbit mAb
Accession(Primary)	O00499
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 65 kDa; Observed MW: 45-80 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:500; ICC/IF-1:100-1:200; FC-1:20-1:50; IP-1:20-1:50

Target information (O00499)

Gene Name	BIN1
About Gene Symbol	BIN1 (Bridging Integrator 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of BIN1 in A549 lysates
using BIN1 antibody.