

BMP2 Rabbit mAb


ABZ2855

Product Overview

Name	BMP2 Rabbit mAb
Accession(Primary)	P12643
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:50-1:100; FC-1:200-1:1000

Target information (P12643)

Gene Name	BMP2
About Gene Symbol	



Western blot analysis of BMP2 in MCF-7, HeLa
lysates using BMP2 antibody.