

BMP4 Rabbit mAb


ABZ2856

Product Overview


Name	BMP4 Rabbit mAb
Accession(Primary)	P12644
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:100-1:500; IP-1:10-1:100

Target information (P12644)

Gene Name	BMP4
About Gene Symbol	



Immunocytochemistry analysis of BMP4 (green) in HEPG2 using BMP4 antibody, and DAPI (blue).



Western blot analysis of BMP4 in K562, HeLa lysates using BMP4 antibody.