


BNIP3L Rabbit mAb

ABZ2857


Product Overview

Name	BNIP3L Rabbit mAb
Accession(Primary)	O60238
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (O60238)

Gene Name BNIP3L**About Gene Symbol**

Western blot analysis of BNIP3L in Jurkat lysates using BNIP3L antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using BNIP3L antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.