

Bora Rabbit mAb


ABZ2858

Product Overview


Name	Bora Rabbit mAb
Accession(Primary)	Q6PGQ7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 61 kDa; Observed MW: 61 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; IP-1:20-1:50

Target information (Q6PGQ7)

Gene Name	BORA
About Gene Symbol	



Western blot analysis of Bora in mouse liver lysates using Bora antibody.



Western blot analysis of Bora in K562, C6 lysates using Bora antibody.