

C4BPB Rabbit mAb


ABZ2860

Product Overview

Name	C4BPB Rabbit mAb
Accession(Primary)	P20851
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 28 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:200; IP-1:50-1:100

Target information (P20851)

Gene Name	C4BPB
About Gene Symbol	C4BPB (Complement Component 4 Binding Protein Beta) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of C4 binding protein in Human serum lysates using C4 Binding Protein antibody.