


C9orf167 Rabbit mAb

ABZ2861

Product Overview

Name	C9orf167 Rabbit mAb
Accession(Primary)	Q9NXH8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (Q9NXH8)

Gene Name TOR4A**About Gene Symbol**

Western blot analysis of C9orf167 in K562, HeLa lysates using C9orf167 antibody.

