

Calbindin Rabbit mAb


ABZ2862

Product Overview


Name	Calbindin Rabbit mAb
Accession(Primary)	P12658
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 30 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; IP-1:20-1:50

Target information (P12658)

Gene Name	Calb1
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded mouse kidney using Calbindin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Calbindin in rat Brain lysates using Calbindin antibody.