

Calretinin Rabbit mAb


ABZ2863

Product Overview


Name	Calretinin Rabbit mAb
Accession(Primary)	P22676
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 32 kDa; Observed MW: 32 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; ICC/IF-1:100-1:200; FC-1:20-1:50

Target information (P22676)

Gene Name	CALB2
About Gene Symbol	



Western blot analysis of Calretinin in rat Brain lysates using Calretinin antibody.



Immunocytochemistry analysis of Calretinin(green) in Hela using Calretinin antibody, and DAPI(blue)