

Carbonic Anhydrase 1 Rabbit mAb


ABZ2864

Product Overview


Name	Carbonic Anhydrase 1 Rabbit mAb
Accession(Primary)	P00915
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 29 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:50-1:100; IHC-1:50-1:100

Target information (P00915)


Gene Name	CA1
About Gene Symbol	



Western blot analysis of Carbonic Anhydrase 1/CA1 in TF1, mouse gastric lysates using Carbonic Anhydrase 1 antibody.



Western blot analysis of Carbonic Anhydrase 1/CA1 in rat Brain lysates using Carbonic Anhydrase 1/CA1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Carbonic Anhydrase I antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.