



CBL Rabbit mAb

ABZ2865


Product Overview

Name	CBL Rabbit mAb
Accession(Primary)	P22681
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 100 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:20-1:50; FC-1:20-1:50

Target information (P22681)

Gene Name CBL**About Gene Symbol**

Immunocytochemistry analysis of CBL (green) in K562 using CBL antibody, and DAPI(blue).



Western blot analysis of CBL in HeLa lysates using CBL antibody.