


CBR3 Rabbit mAb

ABZ2866

Product Overview

Name	CBR3 Rabbit mAb
Accession(Primary)	O75828
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 31 kDa; Observed MW: 31 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:50-1:100

Target information (O75828)

Gene Name CBR3**About Gene Symbol**

Western blot analysis of CBR3 in HeLa lysates using CBR3 antibody.

