

CD11a Rabbit mAb


ABZ2867

Product Overview

Name	CD11a Rabbit mAb
Accession(Primary)	P20701
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 129 kDa; Observed MW: 150 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:50-1:100; FC-1:20-1:50; IP-1:20-1:50

Target information (P20701)

Gene Name	ITGAL
About Gene Symbol	



Western blot analysis of CD11a in THP-1 lysates using Integrin alpha L antibody.