


CD127 Rabbit mAb

ABZ2868

Product Overview

Name	CD127 Rabbit mAb
Accession(Primary)	P16871
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 52 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200

Target information (P16871)

Gene Name IL7R**About Gene Symbol**

Western blot analysis of CD127 in K562 lysates using CD127 antibody.

