

CD22 Rabbit mAb

ABZ2869

Product Overview


Name	CD22 Rabbit mAb
Accession(Primary)	P20273
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 95 kDa; Observed MW: 140 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:50-1:100; IP-1:20-1:50

Target information (P20273)

Gene Name	CD22
About Gene Symbol	



Immunocytochemistry analysis of CD22 (green) in HL-60 using CD22 antibody, and DAPI (blue).



Western blot analysis of CD22 in Raji, K562, Ramos lysates using CD22 antibody.