

## CD31 Rabbit mAb


ABZ2871

### Product Overview


|                             |                                                            |
|-----------------------------|------------------------------------------------------------|
| <b>Name</b>                 | CD31 Rabbit mAb                                            |
| <b>Accession(Primary)</b>   | P16284                                                     |
| <b>Host Species</b>         | Rabbit                                                     |
| <b>Clonality</b>            | Monoclonal                                                 |
| <b>Isotype</b>              | IgG                                                        |
| <b>Species Reactivity</b>   | Human                                                      |
| <b>Calculated MW</b>        | Calculated MW: 83 kDa; Observed MW: 130 kDa                |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:100-1:500;<br>ICC/IF-1:200-1:500 |

### Target information (P16284)

|                          |        |
|--------------------------|--------|
| <b>Gene Name</b>         | PECAM1 |
| <b>About Gene Symbol</b> |        |



Western blot analysis of CD31 in Jurkat lysates using CD31 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using CD31 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.