


CD35 Rabbit mAb

ABZ2873

Product Overview

Name	CD35 Rabbit mAb
Accession(Primary)	P17927
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 224 kDa; Observed MW: 280 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (P17927)

Gene Name CR1**About Gene Symbol**

Western blot analysis of CD35 in Human tonsil lysates using CD35 antibody.

