

CD3D Rabbit mAb


ABZ2874

Product Overview

Name	CD3D Rabbit mAb
Accession(Primary)	P04234
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 19 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:200-1:1000; FC-1:20-1:50; IP-1:50-1:100

Target information (P04234)

Gene Name	CD3D
About Gene Symbol	



Western blot analysis of CD3D in Jurkat lysates
using CD3D antibody.