

## CD55 Rabbit mAb


ABZ2877

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | CD55 Rabbit mAb                            |
| <b>Accession(Primary)</b>   | P08174                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human                                      |
| <b>Calculated MW</b>        | Calculated MW: 41 kDa; Observed MW: 78 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-1:50-1:100            |

### Target information (P08174)

|                          |                                                                                                                                                                                                                                                                                       |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | CD55                                                                                                                                                                                                                                                                                  |
| <b>About Gene Symbol</b> | CD55 (CD55 Molecule (Cromer Blood Group)) is a cell surface marker involved in adaptive immunity. It is an important research target in immunology, flow cytometry, and cancer immunotherapy and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of CD55 in HeLa lysates  
using CD55 antibody.